



Parainfluenza-Virus Typ 3-Antikörper

Testparameter anti-Parainfluenza-Virus 3-IgG (Serum)
anti-Parainfluenza-Virus 3-IgG (Liquor)

Material Serum, 1 mL
Liquor, 1 mL

Referenzbereich	Serum	[U/mL]	Grauzone
	anti-Parainfluenza-Virus 3-IgG	< 65	65 - 95
	Liquor		
	anti-Parainfluenza-Virus 3-IgG	< 65	
	ASI < 2		

Methode Elisa

Qualitätskontrolle intern

Anforderungsschein Download und Analysenposition

Auskünfte Infektionsimmunologie

Analysenkosten EBM, GOÄ