



Rickettsien-Antikörper

Testparameter anti-Rickettsia conorii / rickettsii -IgG
anti-Rickettsia conorii / rickettsii -IgM
anti-Rickettsia mooseri / typhi -IgG
anti-Rickettsia mooseri / typhi -IgM

Material Serum, 1 mL

Referenzbereich	[IIF-TE]
anti-Rickettsia conorii / rickettsii -IgG	IgG: 1 : < 40
anti-Rickettsia conorii / rickettsii -IgM	IgM: 1 : < 20
anti-Rickettsia mooseri / typhi -IgG	IgG: 1 : < 40
anti-Rickettsia mooseri / typhi -IgM	IgM: 1 : < 20

Methode IIFT

Qualitätskontrolle intern

Anforderungsschein Download und Analysenposition

Auskünfte Infektionsimmunologie

Analysenkosten EBM, GOÄ